

Terumo® Syringe with Needle



Product Information Card

Description:

Sterile 3-part syringe for single use pre-connected with a sterile hypodermic needle. Sterile hypodermic syringes are intended for the aspiration of fluid or blood, or for the injection of fluids immediately after filling into the human body.

Components:

Description	Raw Materials
Syringe	
Barrel	Polypropylene
Plunger	Polypropylene
Gasket	Thermoplastic elastomer
Lubricant	Silicone oil
Needle	
Hub	Polypropylene + master batch
Cannula	Stainless steel
Adhesive	Epoxy glue
Lubricant	Silicone



Product Range:

Product Code	Syringe Scale			Cannula				Replacement for
	Nominal	Extended	Increments	Ø	Length			
SS+01T25161	1 ml	1 ml	0.01	25G	0.50 mm	5/8 inch	16 mm	BS-01T2516
SS+01T26131	1 ml	1 ml	0.01	26G	0.45 mm	1/2 inch	12 mm	BS-01T2613
SS-02S2116	2 ml	2.5 ml	0.1	21G	0.80 mm	5/8 inch	16 mm	
SS+02S21251	2 ml	2.5 ml	0.1	21G	0.80 mm	1 inch	25 mm	SS-02S2125
SS+02S21381	2 ml	2.5 ml	0.1	21G	0.80 mm	1 1/2 inch	40 mm	SS-02S2138
SS+02S22381	2 ml	2.5 ml	0.1	22G	0.70 mm	1 1/2 inch	40 mm	SS-02S2238
SS+02S23251	2 ml	2.5 ml	0.1	23G	0.60 mm	1 inch	25 mm	SS-02S2325
SS+02S23321	2 ml	2.5 ml	0.1	23G	0.60 mm	1 1/4 inch	30 mm	SS-02S2332
SS-02S2516	2 ml	2.5 ml	0.1	25G	0.50 mm	5/8 inch	16 mm	
SS+05S21381	5 ml	6 ml	0.2	21G	0.80 mm	1 1/2 inch	40 mm	SS*05SE21381
SS+05S22381	5 ml	6 ml	0.2	22G	0.70 mm	1 1/2 inch	40 mm	SS*05SE22381
SS+05S23321	5 ml	6 ml	0.2	23G	0.60 mm	1 1/4 inch	30 mm	SS*05SE23321
SS+10S20381	10 ml	12 ml	0.2	20G	0.90 mm	1 1/2 inch	40 mm	SS-10S2038
SS+10S21381	10 ml	12 ml	0.2	21G	0.80 mm	1 1/2 inch	40 mm	SS-10S2138
SS+10S22381	10 ml	12 ml	0.2	22G	0.70 mm	1 1/2 inch	40 mm	SS-10S2238

Sterilization:

For - code items: with ethylene oxide to a Sterility Assurance Level of at least 10⁻⁶
 For * and + code items: with e-beam radiation to a Sterility Assurance Level of at least 10⁻⁶

Storage:

Shelf life of 5 years

* The lot number and expiry date are printed on the paper or embossed in the film of the blister pack

* Do not store at extreme temperature and humidity

Terumo® Syringe with Needle

Packaging:

Single unit: blister pack

Product Code	Unit Box	Shipping Carton		
		Quantity	Gross Weight	Dimensions (LxWxH)
SS+01T25161	100	1800	9.8 kg	640 x 366 x 335 mm
SS+01T26131	100	1800	9.8 kg	640 x 366 x 335 mm
SS-02S2116	100	600	3.6 kg	552 x 238 x 302 mm
SS+02S21251	100	1800	9.7 kg	634 x 344 x 377 mm
SS+02S21381	100	1800	9.7 kg	634 x 344 x 377 mm
SS+02S22381	100	1800	9.7 kg	634 x 344 x 377 mm
SS+02S23251	100	1800	9.7 kg	634 x 344 x 377 mm
SS+02S23321	100	1800	9.7 kg	634 x 344 x 377 mm
SS-02S2516	100	600	3.6 kg	552 x 238 x 302 mm
SS+05S21381	100	1200	9.4 kg	480 x 364 x 485 mm
SS+05S22381	100	1200	9.4 kg	480 x 364 x 485 mm
SS+05S23321	100	1200	9.4 kg	480 x 364 x 485 mm
SS+10S20381	100	800	9.6 kg	540 x 408 x 365 mm
SS+10S21381	100	800	9.6 kg	540 x 408 x 365 mm
SS+10S22381	100	800	9.6 kg	540 x 408 x 365 mm

Precautions for Use: Check integrity of the unit package before use. For single use only. Dispose of safely after single use to avoid risk of infection.

Status of Regulation: Class IIa, medical device following MDD 93/42/EEC

Notified body: TÜV Rheinland 0197

For * code items:

Manufacturer: TERUMO CORPORATION JAPAN
EC representative: TERUMO EUROPE NV

For + code items:

Manufacturer: TERUMO (PHILIPPINES) CORPORATION
EC representative: TERUMO EUROPE NV

For - code items:

Manufacturer: TERUMO EUROPE NV
Interleuvenlaan 40
B - 3001 Leuven
Belgium



Terumo Corporation
44-1, 2-chome, Hatagaya
Shibuya-ku, Tokyo 151-0072,
Japan
Tel.: +81 3 3374 8111
Fax: +81 3 3374 8196

Terumo Europe NV
Interleuvenlaan 40
3001 Leuven, Belgium
Tel.: +32 16 38 12 11
Fax: +32 16 40 02 49

Terumo Europe NV
Benelux Sales Division
Belgium:
Tel.: +32 16 39 25 80
The Netherlands:
Tel.: 0800 0220396

Terumo Europe NV
Emerging Markets Division
Tel.: +32 16 38 13 08

Terumo Deutschland GmbH
Tel.: +49 6196 80 230

Terumo Deutschland GmbH
Zweigniederlassung Spreitenbach
Tel.: +41 56 419 10 10

Laboratoires Terumo France S.A.S.
Tel.: +33 130 96 13 00

Terumo Italia S.r.l.
Tel.: +39 0651 96 14 20

Terumo Europe España SL
Tel.: +34 902 10 12 98

Terumo UK Ltd
Tel.: +44 1276 480 440

Terumo Russia LLC
Tel.: + 7 495 988 4740

Terumo Sweden AB
Tel.: +46 3174 85 880

Terumo Denmark
Tel.: +45 7020 93 80

Terumo Corporation Dubai Branch
(Middle East and Africa)
Tel.: + 971 4 221 2220



HO 76GB-0908TE-VII(12.14)E

© Registered Trademark
Published by Terumo Europe N.V.